



Summary validation Delvotest® SP NT and Delvotest® T by Laboratories MUVA Kempten GmbH (Germany) and ARA Sardegna (Italy)

DSM Food Specialties B.V.

Alexander Fleminglaan 1
2613 AX Delft
P.O. Box 1
2600 MA Delft
Netherlands

phone +31(0)297 624619
fax +31(0)297 624619

Date
22 March, 2017

TO WHOM IT MAY CONCERN

DSM commissioned an independent study executed by the following laboratories: MUVA Kempten and ARA Sardegna (Associazione Regionale Allevatori della Sardegna); which are ISO/IEC 17025 certified, with the aim to verify the limits of detection (LOD) of Delvotest® SP NT and Delvotest® T in both ampoules and plates variants.

Delvotest® SP NT and Delvotest® T tests were used as qualitative screening methods for detecting antibiotic residues in milk.

The test matrix was raw bovine milk with incubation at control time (CT). The test end-point was measured both visually and instrumentally with Delvo®Scan.

The sensitivity was tested by using blank and spiked samples at levels of interest.

As per requested, below is the sensitivity of Delvotest® SP NT and Delvotest® T, found by above mentioned laboratories, with respect to the following substances:

Although diligent care has been used to ensure that the information provided herein is accurate, nothing contained herein can be construed to imply any representation or warranty for which we assume legal responsibility, including without limitation any warranties as to the accuracy, currency or completeness of this information or of non-infringement of third party intellectual property rights. The content of this document is subject to change without further notice. This document is non-controlled and will not be automatically replaced when changed. Please contact us for the latest version of this document or for further information. Since the user's product formulations, specific use applications and conditions of use are beyond our control, we make no warranty or representation regarding the results which may be obtained by the user. It shall be the responsibility of the user to determine the suitability of our products for the user's specific purposes and the legal status for the user's intended use of our products.

DSM Food Specialties B.V. | A. Fleminglaan 1 | 2613 AX Delft | The Netherlands | info.food@dsm.com | Trade Register Number 27235314



Delvotest® SP NT	Sensitivity (ppb)		
	Antibiotic	MRL	Plates
Sulfadiazine	100	65	40
Meloxicam	15	*Neg. at 1500	*Neg. at 1500
Trimethoprim	50	160	160
Albendazole	NA	*Neg. at 5000	*Neg. at 5000

Delvotest® T	Sensitivity (ppb)		
	Antibiotic	MRL	Plates
Sulfadiazine	100	40	40
Meloxicam	15	*Neg. at 1500	*Neg. at 1500
Trimethoprim	50	100	100
Albendazole	NA	*Neg. at 5000	*Neg. at 5000

**Substances not tested beyond the mentioned value.*

on behalf of DSM Food Specialties B.V.

Bibi Duyvesteyn
Business Line Manager - Delvotest®

Although diligent care has been used to ensure that the information provided herein is accurate, nothing contained herein can be construed to imply any representation or warranty for which we assume legal responsibility, including without limitation any warranties as to the accuracy, currency or completeness of this information or of non-infringement of third party intellectual property rights. The content of this document is subject to change without further notice. This document is non-controlled and will not be automatically replaced when changed. Please contact us for the latest version of this document or for further information. Since the user's product formulations, specific use applications and conditions of use are beyond our control, we make no warranty or representation regarding the results which may be obtained by the user. It shall be the responsibility of the user to determine the suitability of our products for the user's specific purposes and the legal status for the user's intended use of our products.

DSM Food Specialties B.V. | A. Fleminglaan 1 | 2613 AX Delft | The Netherlands | info.food@dsm.com | Trade Register Number 27235314